# **An Endothermic Reaction**

## Purpose

To demonstrate an endothermic reaction by combining two solids.

## Materials

| 250 mL beaker | hydrated barium hydroxide |
|---------------|---------------------------|
| wooden splint | ammonium thiocyanate      |
| wooden board  |                           |

## Procedure

- 1. Place about 20 grams of Ba(OH)<sub>2</sub> 8 H<sub>2</sub>O crystals in a 250 mL beaker.
- 2. Add 10 grams of ammonium thiocyanate to the beaker.
- 3. Stir the two solids together with a wooden splint.
- 4. Place a small pool of water on a piece of board. Place the beaker on the pool of water.
- 5. After 2-3 minutes, the beaker will freeze to the board.

### **Additional Information**

- 1. Ammonium chloride (7.0 grams) or ammonium nitrate (10 grams) may replace the ammonium thiocyanate.
- 2. Ba(OH)<sub>2</sub> 8 H<sub>2</sub>O + 2 NH<sub>4</sub>SCN Ba(SCN)<sub>2</sub> + 2 NH<sub>3</sub> + 10 H<sub>2</sub>O
- A variation involves adding about 10-15 grams of ammonium nitrate to 100 mL of water. The temperature will drop significantly.

### Disposal

Excess solutions/solids should be placed in a properly labeled waste container with UI# 205437.

### Reference

Summerlin, L., Ealy, J.; Chemical Demonstrations: A Sourcebook for Teachers, 1985.

Reactions: Endothermic 25